

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



Thioacetazone

Item No. 34589

CAS Registry No.:	104-06-3	
Formal Name:	N-[4-[[2-(aminothioxomethyl)	н
	hydrazinylidene]methyl]phenyl]-acetamide	
Synonyms:	Thiacetazone, NSC 3550, SQ 2321	N N
MF:	C ₁₀ H ₁₂ N ₄ OS	
FW:	236.3	
Purity:	≥98%	· · · · · · · · · · · · · · · · · · ·
UV/Vis.:	λ _{max} : 330 nm	Η̈́
Supplied as:	A solid	
Storage:	-20°C	
Stability:	≥2 years	
Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.		

Laboratory Procedures

Thioacetazone is supplied as a solid. A stock solution may be made by dissolving the thioacetazone in the solvent of choice, which should be purged with an inert gas. Thioacetazone is soluble in organic solvents such as DMSO and dimethyl formamide. The solubility of thioacetazone in these solvents is approximately 30 mg/ml.

Thioacetazone is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, thioacetazone should first be dissolved in DMSO and then diluted with the aqueous buffer of choice. Thioacetazone has a solubility of approximately 0.25 mg/ml in a 1:3 solution of DMSO:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

Thioacetazone is an antitubercular agent.¹ It is active against M. tuberculosis H37Rv, M. kansasii Bostrom, and M. avium serotype 2 (MICs = 0.3, 0.6, and 20 µg/ml, respectively) but not M. simiae N29 or M. intracellulare serotype 16 (MICs = >40 μ g/ml for both). Formulations containing thioacetazone were previously used in the treatment of tuberculosis.

Reference

1. Collins, F.M., Klayman, D.L., and Morrison, N.E. Activity of 2-acetylpyridine and 2-acetylquinoline thiosemicarbazones tested in vitro in combination with other antituberculous drugs. Am. Rev. Respir. Dis. 125(1), 58-60 (1982).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 08/24/2021

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM