



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



# PRODUCT INFORMATION

## Rp-Adenosine-5'-O-(1-thiotriphosphate) (sodium salt)

Item No. 38535

**Formal Name:** Adenosine, 5'→P''-ester with [P''(R)]-thiotriphosphoric acid ((HO)2P(O)OP(O)(OH)OP(O)(OH)(SH)), tetrasodium salt

**Synonym:** Rp-ATP-α-S

**MF:** C<sub>10</sub>H<sub>12</sub>N<sub>5</sub>O<sub>12</sub>P<sub>3</sub>S • 4Na


**FW:** 611.2

**Purity:** ≥95%

**Supplied as:** A solution in water

**Storage:** -80°C

**Stability:** ≥2 years



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

### Description

Rp-ATP-α-S is an isomer of the sulfur-containing nucleotide derivative ATP-α-S and an agonist of the purinergic P2Y<sub>1</sub> receptor.<sup>1</sup> It increases calcium mobilization in HEK293 cells expressing the human P2Y<sub>1</sub> receptor (EC<sub>50</sub> = 75 nM). Rp-ATP-α-S binds to washed isolated human platelets (K<sub>i</sub> = 156 nM), inhibits ADP-induced aggregation of isolated human platelet-rich plasma (PRP; pA<sub>2</sub> = 4.74), and inhibits cAMP production induced by prostaglandin E<sub>1</sub> (PGE<sub>1</sub>; Item No. 13010) in isolated human PRP (pA<sub>2</sub> = 5.26).<sup>1,2</sup> It also induces relaxation of carbamoylcholine-precontracted guinea pig taenia coli strips (EC<sub>50</sub> = 56 nM).<sup>3</sup> Rp-ATP-α-S has been used in the synthesis of cyclic dinucleotides that are recognized by bacterial riboswitches.<sup>4</sup>

### References

1. Major, D.T., Nahum, V., Wang, Y., *et al.* Molecular recognition in purinergic receptors. 2. Diastereoselectivity of the h-P2Y<sub>1</sub>-receptor. *J. Med. Chem.* **47**(18), 4405-4416 (2004).
2. Agarwal, A.K., Tandon, N.N., Greco, N.J., *et al.* Evaluation of the binding to fixed platelets of agonists and antagonists of ADP-induced aggregation. *Thromb. Haemost.* **62**(4), 1103-1106 (1989).
3. Cusack, N.J. and Hourani, S.M. Adenosine 5-diphosphate antagonists and human platelets: No evidence that aggregation and inhibition of stimulated adenylate cyclase are mediated by different receptors. *Br. J. Pharmacol.* **76**(1), 221-227 (1982).
4. Burnstock, G., Cusack, N.J., and Meldrum, L.A. Effects of phosphorothioate analogues of ATP, ADP and AMP on guinea-pig taenia coli and urinary bladder. *Br. J. Pharmacol.* **82**(2), 369-374 (1984).
5. Launer-Felty, K.D. and Strobel, S.A. Enzymatic synthesis of cyclic dinucleotide analogs by a promiscuous cyclic-AMP-GMP synthetase and analysis of cyclic dinucleotide responsive riboswitches. *Nucleic Acids Res.* **46**(6), 2765-2776 (2018).

#### WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

#### SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the [complete](#) Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 06/27/2023

#### CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD  
ANN ARBOR, MI 48108 · USA

**PHONE:** [800] 364-9897  
[734] 971-3335

**FAX:** [734] 971-3640

CUSTSERV@CAYMANCHEM.COM  
[WWW.CAYMANCHEM.COM](http://WWW.CAYMANCHEM.COM)