

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



Ferritin PROTAC DeFer-2

Item No. 38931

Formal Name:	(2S,4R)-1-((S)-3,3-dimethyl-2-(4- oleamidobutanamido)butanoyl)- 4-hydroxy-N-(4-(4-methylthiazol- 5-yl)benzyl)pyrrolidine-2- carboxamide	
Synonym:	Ferritin Proteolysis-targeting	H O N
	Chimera DeFer2	
MF:	C ₄₄ H ₆₉ N ₅ O ₅ S	
FW:	780.1	
Purity:	≥98%	
Supplied as:	A solid	
Storage:	-20°C	
Stability:	≥3 years	
Information represents the product exections. Databased in a line republic are provided on each continues of an above		

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Ferritin PROTAC DeFer-2 is supplied as a solid. A stock solution may be made by dissolving the ferritin PROTAC DeFer-2 in the solvent of choice, which should be purged with an inert gas. Ferritin PROTAC DeFer-2 is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide. The solubility of ferritin PROTAC DeFer-2 in these solvents is approximately 30, 10, and 20 mg/ml, respectively.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of ferritin PROTAC DeFer-2 can be prepared by directly dissolving the solid in aqueous buffers. Ferritin PROTAC DeFer-2 is slightly soluble in PBS (pH 7.2). We do not recommend storing the aqueous solution for more than one day.

Description

Ferritin PROTAC DeFer-2 is a proteolysis-targeting chimera (PROTAC) containing oleic acid (Item Nos. 90260 | 24659), which binds to ferritin dimers, conjugated to von Hippel-Lindau ligand (VHL ligand 1; Item No. 21591), which binds to VHL E3 ligase, via an alkyl linker.¹ It induces ferritin degradation in B16/F10 murine melanoma cells when used at a concentration of 0.5 μ M, an effect not induced by a negative control with an inverted VHL ligand 1 stereocenter, indicating that the E3 ligase activity of VHL ligand 1 is required for this effect. Ferritin PROTAC DeFer-2 (5 μM) increases intracellular levels of free iron and induces pyroptosis in B16/F10 cells. Ferritin PROTAC De-Fer-2 (10 mg/kg) encapsulated in albumin nanoparticles reduces tumor growth and increases survival in a B16/F10 murine melanoma model.

Reference

1. Chen, Y., Li, W., Kwon, S., et al. Small-molecule ferritin degrader as a pyroptosis inducer. J. Am. Chem. Soc. 145(17), 9815-9824 (2023).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 07/19/2023

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM