

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



C6 NBD Ceramide (d18:1/6:0)

Item No. 62527

CAS Registry No.:	94885-02-6		
Formal Name:	N-[7-(4-nitrobenzo-2-oxa-1,3-diazole)]-		
	6-aminocaproyl-D-erythro-sphingosine	NO ₂	
Synonyms:	C6 NBD Ceramide, NBD Cer(18:1/2:0),	N	Ļ
	NBD Ceramide (d18:1/6:0), N-C6-		$\hat{\mathbf{o}}$
	NBD-D-erythro-Sphingosine		ŕ
MF:	$C_{30}H_{49}N_5O_6$	0	
FW:	575.7		
Purity:	≥98%		
UV/Vis.:	λ _{max} : 229, 333, 465 nm	ОН	
Supplied as:	A crystalline solid	ОН	
Storage:	-20°C		
Stability:	≥4 years		

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

C6 NBD ceramide (d18:1/6:0) is supplied as a crystalline solid. A stock solution may be made by dissolving the C6 NBD ceramide (d18:1/6:0) in the solvent of choice. C6 NBD ceramide (d18:1/6:0) is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF), which should be purged with an inert gas. The solubility of C6 NBD ceramide (d18:1/6:0) in these solvents is approximately 5, 10, and 25 mg/ml, respectively.

C6 NBD ceramide (d18:1/6:0) is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, C6 NBD ceramide (d18:1/6:0) should first be dissolved in DMF and then diluted with the aqueous buffer of choice. C6 NBD ceramide (d18:1/6:0) has a solubility of approximately 500 μ g/ml in a 1:1 solution of DMF:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

C6 NBD ceramide is a biologically active fluorescent analog of short-chain, membrane-permeable ceramides.¹ It is as effective as C6 ceramide in the inhibition of of viral glycoprotein transport through the Golgi.¹ C6 NBD ceramide has been used as a fluorescent substrate for the activity of UDP-glucose:ceramide glycosyltransferase and to demonstrate the translocation of glucocerebroside and sphingomyelin from the Golgi to the plasma membrane.^{2,3}

References

- 1. Rosenwald, A.G. and Pagano, R.E. Inhibition of glycoprotein traffic through the secretory pathway by ceramide. J. Biol. Chem. 268, 4577-4579 (1993).
- 2. Paul, P., Kamisaka, Y., Marks, D.L., et al. Purification and characterization of UDP-glucose: Ceramide glucosyltransferase from rat liver Golgi membranes. J. Biol. Chem. 271, 2287-2293 (1996).
- 3. Lipsky, N.G. and Pagano, R.E. Intracellular translocation of fluorescent sphingolipids in cultured fibroblasts: Endogenously synthesized sphingomyelin and glucocerebroside analogues pass through the golgi apparatus en route to the plasma membrane. J. Cell. Biol. 100, 27-34 (1985).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 03/28/2024

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM