

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Linoleic Acid-d₁₁ Item No. 9002193

Formal Name: 9Z.12Z-octadecadieno

ic-14,14',15,15',16,16',17,17',18,18,18-

Synonyms: Telfairic Acid-d₁₁, C18:2(9Z,12Z)-d₁₁,

C18:2 n-6-d11, 9, 12-Octadecadienoic

Acid-d₁₁

 $C_{18}H_{21}\bar{D}_{11}O_2$ MF:

291.5 FW:

Chemical Purity: ≥98% (Linoleic Acid)

Deuterium

Incorporation: \geq 99% deuterated forms (d₁-d₁₁); \leq 1% d₀

Supplied as: A solution in methyl acetate

Storage: -20°C Stability: ≥1 year

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



 $\ \ \text{Linoleic acid-d}_{11} \ \text{is supplied as a solution in methyl acetate. To change the solvent, simply evaporate the } \\$ linoleic acid- d_{11} under a gentle stream of nitrogen and immediately add the solvent of choice. Solvents such as ethanol, DMSO, and dimethyl formamide purged with an inert gas can be used. The solubility of linoleic acid-d₁₁ in these solvents is at least 100 mg/ml.

Linoleic acid-d₁₁ is used as an internal standard for the quantification of linoleic acid by stable isotope dilution MS. The accuracy of the sample weight in this vial is between 5% over and 2% under the amount shown on the vial. If better precision is required, the deuterated standard should be quantitated against a more precisely weighed unlabeled standard by constructing a standard curve of peak intensity ratios (deuterated versus unlabeled).

Description

Linoleic acid-d₁₁ is intended for use as an internal standard for the quantification of linoleic acid (Item Nos. $90150 \mid 90150.1 \mid 21909$) by GC- or LC-MS. Linoleic acid is an essential ω -6 polyunsaturated fatty acid (PUFA).¹ It is the most abundant PUFA in a variety of foods, and dietary sources of linoleic acid include vegetable oils, meats, nuts, seeds, and eggs. Linoleic acid (30 μM) increases migration of IEC-6 rat intestinal epithelial cells in a wound healing assay. 2 Rats fed a linoleate-deficient diet exhibit decreased body weight and an increased ratio of eicosatrienoate to arachidonate in liver and serum phospholipids compared with rats fed a control diet, as well as mild scaling of forepaw skin.³

References

- 1. Whelan, J., Fritsche, K. Linoleic acid. Adv. Nutr. 4(3), 311-312 (2013).
- Ruthig, D.J., Meckling-Gill, K.A. Both (n-3) and (n-6) fatty acids stimulate wound healing in the rat intestinal epithelial cell line, IEC-6. J. Nutr. 129(10), 1791-1798 (1999).
- 3. Cunnane, S.C., Anderson, M.J. Pure linoleate deficiency in the rat: Influence on growth, accumulation of n-6 polyunsaturates, and [1-14C]linoleate oxidation. J. Lipid Res. 38(4), 805-812 (1997).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 09/10/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM