



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

POOLED NORMAL RABBIT SERUM

CL1000
CL1000-10
CL1000-20
CL1000-100
CL1000-500/ 500S
CL1000-1000

DESCRIPTION:

Cedarlane Laboratories has developed a normal rabbit serum specifically for use in immunohistopathology. This serum has been processed to minimize haemolysis and filtered to 0.45µm (non-sterile), 0.22µm (sterile) to remove protein aggregates. CL1000 is preserved with 0.1% sodium azide to allow for long term storage at 4°C.

PRESENTATION:

CL1000: 10 ml, liquid with 0.1% NaN₃, filtered 0.22µm sterile.
CL1000-10: 10ml, frozen, filtered 0.45µm, non-sterile.
CL1000-20: 20 ml, frozen, filtered 0.45µm, non-sterile.
CL1000-100: 100ml, frozen, filtered 0.45µm, non-sterile.
CL1000-500: 500ml, frozen, filtered 0.45µm, non-sterile.
CL1000-500S: 500ml, frozen, filtered 0.22µm, sterile.
CL1000-1000: 1000ml, frozen, filtered 0.45µm, non-sterile.

STORAGE/STABILITY:

Store CL1000 at 4°C.

Note: Sera may precipitate in storage.

Store CL1000-10; CL1000-20; CL1000-100; CL1000-500; CL1000-500S; CL1000-1000 at -20°C or below.

APPLICATIONS:

Blocking sera, control sera, a serum additive in dilution buffers.

FOR RESEARCH USE ONLY

JCr 01/25/06

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com