



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## NORMAL SWINE SERUM

**CL1100**  
**CL1100-30**  
**CL1100-100**

### **DESCRIPTION:**

Cedarlane Laboratories has developed a normal swine serum specifically for use in immunohistopathology. It has been processed for minimum haemolysis and filtered to 0.22µm. CL1100 is preserved with 0.1% sodium azide to allow for long term storage at 4°C.

### **PRESENTATION:**

**CL1100:** 10 ml, liquid with 0.1% NaN<sub>3</sub>.  
**CL1100-30:** 30ml, frozen, filtered 0.22µm, sterile.  
**CL1100-100:** 100ml, frozen, filtered 0.22µm, sterile.

### **STORAGE/STABILITY:**

Store CL1100 at 4°C.  
**Note: Sera may precipitate in storage.**  
Store CL1100-30 and CL1100-100 at -20°C or below.

### **APPLICATIONS:**

Blocking sera, control sera, a serum additive in dilution buffers.

**FOR RESEARCH USE ONLY**

JK 032812

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)