



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

TECHNICALLY Speaking

Place your order with CEDARLANE® or your local distributor.
Please contact CEDARLANE® for lot specific information.

Sheep Anti-Rat von Willebrand Factor Purified IgG

CL20176A-R
Lot : AB59-12R2
Exp: 06/2011

- Immunogen:** Purified von Willebrand Factor (vWF) from rat plasma.
- Presentation:** 1.0 ml of a purified preparation of immunoglobulin prepared from sheep antiserum. The IgG concentration is 5 mg/ml in 0.01 M HEPES, pH 7.4 containing 0.15 M NaCl and 50% glycerol (v/v).
- Storage:** Store in the freezer at -20°C. Avoid storage in frost-free freezers. Allow product to warm to room temperature and gently mix before use.
- Specificity:** This antibody is specific for rat vWF as determined by immunoelectrophoresis and ELISA.
- Species Cross Reactivity:** (immunodiffusion vs. citrated plasma)
- | | | | | | |
|---------|----|--------|----|------|----|
| Human: | ND | Mouse: | ND | Rat: | ND |
| Rabbit: | ND | Pig: | ND | Dog: | ND |
- Application:** This product can be used for immunoelectrophoresis and direct ELISA. This product is **not** recommended for use with mouse plasma in ELISA.
- Reference:**
1. Montgomery RR, Collier BS, in Hemostasis and Thrombosis, 3rd Edition, eds. RW Colman, J Hirsh, VJ Marder and EW Salzman, pp. 134-168, J.B. Lippincott Co., Philadelphia PA, USA, 1994.
 2. Sadler JE; von Willebrand Factor (Minireview); JBC 266:34, pp 22777-22780, 1991

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory
Laboratory Reagent For Research Use Only

JV 09/23/08

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com