



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Place your order with CEDARLANE® or your local distributor.  
*Please contact CEDARLANE® for lot specific information.*

## Purified Chicken Anti-Human Hemoglobin

**CL2114AP**

**Lot:**

<b>Description:</b>	Anti-human Hemoglobin was affinity purified on human Hemoglobin agarose.
<b>Antigen:</b>	Human Hemoglobin
<b>Host Animal:</b>	Chicken WL-36
<b>Form:</b>	Liquid. Chicken IgG fraction in 0.15 M NaCl, 0.02 M NaH <sub>2</sub> PO <sub>4</sub> , 0.1% NaN <sub>3</sub> , pH 7.2.
<b>Volume:</b>	2mL
<b>Storage:</b>	4°C - 8°C
<b>Performance:</b>	This antibody can be used in electroblotting with a suggested dilution of 1:1000. This antibody is specific against human hemoglobin in electroblotting when tested against human red blood cell lysate.

**FOR RESEARCH USE ONLY**

TN 11/25/04

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



*toll free: 1-800-268-5058*  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)