



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

STANDARD GUINEA PIG COMPLEMENT

(Pooled Normal Guinea Pig Serum)

**CL5000F-R
CL5000-50-R**

CL5000-100-R

DESCRIPTION:

The serum is intended for use in systems requiring an active source of Guinea Pig complement including hemolytic assays, complement fixation, and hemolytic antibody titrations. Each batch is tested for high activity and low background hemolysis in a standard hemolytic assay.

<u>PRESENTATION:</u>	CL5000F-R	5 ml, frozen
	CL5000-50-R	50 ml, frozen
	CL5000-100-R	100ml, frozen

STORAGE: Store at -70°C

USE:

Keep on ice, and use within 3 hours after thawing. If smaller volumes are required, aliquot to appropriate volumes immediately and store frozen material at -70°C. Upon the second thaw, discard all unused complement. Exposure of complement to more than one freeze/thaw cycle will damage the product.

STERILITY:

This reagent is not sold as sterile, but can be sterilized by filtration if necessary. Dilute to working dilution and then filter to avoid loss of product.

FOR RESEARCH USE ONLY

JK 07/15/19

Visit our website for your local distributor.

CEDARLANE® 
www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com