



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically  
Speaking*

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Llama Complement

**CL6001-10**  
**CL6001-50**

**CL6001-100**  
**CL6001-500**

### **DESCRIPTION:**

Cedarlane's llama complement is freshly prepared from serum. It is processed from whole blood at 0-4°C quickly in order to maintain complement activity. It is intended for use as a source of complement for any of the standard lymphocytotoxicity and hemolytic procedures.

### **PRESENTATION:**

**CL6001-10:** 10 ml, frozen  
**CL6001-50:** 50 ml, frozen  
**CL6000-100:** 100 ml, frozen  
**CL6001-500:** 500 ml, frozen

**STORAGE:** Store at -70°C.

### **STERILITY:**

This reagent is not sold as sterile but can be sterilized by filtration if necessary. Dilute to working dilution and then filter to avoid loss of product.

### **APPLICATIONS:**

This reagent is suitable for use in lymphocytotoxicity and hemolytic procedures.

**FOR RESEARCH USE ONLY**

JK 08/04/20

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2015 and ISO 13485:2016  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)