



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Rabbit anti-Human Transferrin

CLA2240

Lot: 0174B

DESCRIPTION:

Cedarlane's **CLA2240** is rabbit antiserum directed against human transferrin. The antiserum of individual host animals was monitored for titer and specificity before pooling into uniform high titer lots.

PRESENTATION:

2 ml, liquid.

STORAGE:

Divide into useful aliquots and store at -20°C.

PRESERVATIVES:

0.05% sodium azide

APPLICATIONS:

This antibody is suitable for use in immunoprecipitation and immunodiffusion.

OUCHTERLONY: RESULTS

	Test #1
Sample Size:	20 µl neat antisera
Distance Between Wells:	4 mm
Test time:	72 hours at 37°C
Sample tested:	10 λ normal human serum diluted 1:2 to 1:512 in PBS 10 λ human transferrin, purified
Results:	Single homogenous band at all dilutions of normal human serum continuous with transferrin

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Laboratory Reagent For Research Use Only

MW 03/29/16

Visit our website for your local distributor.

CEDARLANE® 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com