



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Anti-Xyloglucan, Monoclonal Antibody

CLAS18-4203-1ML

Immunogen: Neoglycoprotein, Xyloglucan

Host: Rat

Clone: LM15

Isotype: IgG2c

Antibody format: Cell culture supernatant, with 0.05% sodium azide.

Quantity: 1 ml

Storage: Store at +4°C (short term) and at -20°C (long term).

Additional Info: Reacts with the XXXG motif of land plants xyloglucan.
This antibody is directed toward the nonreducing end, with the requirement for the second to last Glc to be substituted with an α -1,6-linked Xyl Ruprecht et al. (2017).

Application Information:

Recommended Dilution: 1:10 (ELISA, IF)

Confirmed Reactivity: Higher plants, ferns and mosses

Not reactive in: No confirmed exceptions from predicted reactivity are currently known.

Selected references: Ruprecht et al. (2017). A Synthetic Glycan Microarray Enables Epitope Mapping of Plant Cell Wall Glycan-Directed Antibodies. *Plant Physiol.* 2017 Nov;175(3):1094-1104. doi: 10.1104/pp.17.00737.

Laboratory Reagent For Research Use Only

JK 06/03/22

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com