



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Affinity Purified HRP Goat anti-Rabbit IgG (H+L)

CLGRHRP

Lot: 40208

Product: Goat Anti-Rabbit (H&L) HRP, Affinity Purified.

Size: 0.5 mg in 0.5 ml buffer containing 0.01 M sodium phosphate, 0.25 M NaCl, pH 7.6 containing 25% glycerol.

Stabilizer: 15 mg/ml BSA

Suggested Dilution Range: 1:250-1:3,000 for IHC/ICC
1:2,000-1:5,000 for ELISA, Western Blot, color substrates
1:5,000-1:10,000 for Western Blot, chemiluminescence

Warning: do not use sodium azide in the conjugate dilution buffer since azide will substantially inhibit horseradish peroxidase (HRP) activity.

Preservative: None.

Short Term Storage: 2-8°C as an undiluted liquid for several weeks. After dilution, do not use for more than one day.

Long Term Storage: -20°C for up to one year.

Specificity: The antibody reacts with both heavy and light chains of rabbit IgG, slight cross-reactivity with IgG of other mammalian species.

Laboratory Reagent For Research Use Only.

JV 09/08/08

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com