



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Superoxide Dismutase Human Recombinant

CLPRO286
CLPRO286-2
CLPRO286-3

Introduction:

Human Cu/Zn Superoxide Dismutase (SOD1) catalyzes the reaction between superoxide anions and hydrogen to yield molecular oxygen and hydrogen peroxide. The enzyme protects the cell against dangerous levels of superoxide. SOD1 binds copper and zinc ions and is 1 of 3 isozymes accountable for destroying free superoxide radicals in the body. The encoded protein neutralizes supercharged oxygen molecules, which can damage cells if their levels are not controlled. Mutations in SOD1 cause a form of familial amyotrophic lateral sclerosis.

Description:

Recombinant Human Cu/Zn Superoxide Dismutase produced in E.Coli is a homodimer, non-glycosylated polypeptide chain containing 2 x 154 amino acids and having a total molecular mass of 31,608 Dalton. The Cu/Zn SOD is purified by proprietary chromatographic techniques.

Synonyms: Superoxide dismutase [Cu-Zn], EC 1.15.1.1, SOD1, SOD, ALS, ALS1, IPOA.

Source: *Escherichia coli*.

Presentation:

20 µg (**CLPRO285**), 100 µg (**CLPRO285-2**), or 1 mg (**CLPRO285-3**) sterile filtered white lyophilized (freeze-dried) powder. Lyophilized from a 0.2 µm filtered concentrated (1.0 mg/ml) solution PBS, pH 7.4

Solubility:

It is recommended to reconstitute the lyophilized SOD in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Stability:

Lyophilized SOD although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution SOD should be stored at 4°C between 2-7 days and for future use below -18°C.

For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).

Please avoid freeze-thaw cycles.

Continued Overleaf.....

Visit our website for your local distributor.

CEDARLANE®
www.cedarlanelabs.com



An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com

Purity:

Greater than 95.0% as determined by:

- a) Analysis by RP-HPLC
- b) Analysis by SDS-PAGE
- c) **Amino Acid Sequence:** The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Thr-Lys-Ala-Val.
- d) **Biological Activity:** The potency per mg was tested by Pyrogalllic Acid method and was found to be more than 10,000 Units/mg.

Laboratory Reagent For Research Use Only