



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Human Serum Albumin (HSA)

CLPRO354
CLPRO354-2
Lot:

Description: Contains 584 amino acid residues derived from the prototypical human serum albumin sequence. Suitable for use in biochemical, excipient (an inert substance used as a diluent or vehicle for a drug), culture media and chromatographic applications.

Source: Human Serum

Presentation: 100 mg (CLPRO354) or 500 mg (CLPRO354-2) sterile filtered clear orange solution at a concentration of 250 mg/ml solution and containing no additives.

Stability: HSA although stable at room temperature for 2 weeks, should be stored at 4°C.

Purity:

Greater than 99.0% as determined by both:

- (a) Analysis by RP-HPLC.
- (b) Analysis by SDS-PAGE.

Applications:

Formulation of Protein Therapeutics
Cell Storage: Cryopreservation
Vaccine formulation and manufacturing
Development of mammalian cell cultures
Infertility treatments
Coating for medical devices
Drug delivery
In vivo diagnostics

Laboratory Reagent for Research Use Only

BA 09/19/12

Visit our website for your local distributor.

CEDARLANE® 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com