



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.
Please contact CEDARLANE® for lot specific information.

Recombinant Human S100A1 Protein

CLPRO364
CLPRO364-2
Lot: 905PS100A101

DESCRIPTION:

S100 belongs to the family of calcium binding proteins with EF-hand type Ca^{+2} binding motive. S100A1 (Calcium Binding Protein A1) is involved in the activation of sarcoplasmic calcium release and the regulation of intermediate filament polymerization.

SOURCE:

Recombinant full length human S100A1 protein expressed in *E.coli*

PRESENTATION:

5 µg (**CLPRO364**) or 10 µg (**CLPRO364-2**) of purified human S100A1 protein at 0.2 mg/ml in 20mM HEPES-KOH, pH7, 50mM NaCl, 1mM EDTA and 1mM DTT.

**APPLICATION AND
SUGGESTED DILUTIONS:**

- ELISA
- Inhibition Assays
- Western Blotting

CHARACTERIZATION:

S100A1 antibody reactive.

ACTIVITY:

Assayed for calcium binding.

PURITY:

Greater than 95% by SDS-PAGE.

STORAGE & STABILITY:

Store vial at -20°C to -80°C . When stored at the recommended temperature, this protein is stable for 12 months.

FOR RESEARCH USE ONLY

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com