



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Recombinant Human S100A11

CLPRO385
CLPRO385-2
Lot: 206S100A11

Description: The S100 family of proteins contain two EF-hand calcium-binding motifs and are thought to be involved in the regulation of a number of cellular processes including cell cycle progression and differentiation. This protein may function in motility, invasion and tubulin polymerisation.
Recombinant Human S100A11 also called S100C is expressed in *E. coli* having a molecular weight of 17kDa fused to an amino terminal hexahistidine tag.

Presentation: 2µg (CLPRO385) or 10µg (CLPRO385-2) rhS100A11 is supplied at a concentration of 1.96 mg/ml in PBS containing 50% glycerol.

Source: *Escherichia Coli*.

Application: rHuS100A11 can be used directly as a positive control in Western blotting, ELISA, immunoprecipitation and other immunological experiments.
The biological activity of this product has not yet been tested.

Storage Conditions: Store at 4°C if entire vial will be used within 1-2 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

Purity:

Greater than 95% (SDS-PAGE).

2 bands on Western blot at 17 and 34 kDa, respectively representing monomeric and dimeric form.

FOR RESEARCH USE ONLY

JV 11/03/06

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com