



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

S100A6

CLPRO409

CLPRO409-2

Lot: 106S100A601

Description: S100A6 also called Calcyclin has been purified by HPLC (see Kuznicki et al. (1989) Biochem. J. 263: 951-956).

Presentation: 2 µg (CLPRO409) or 5µg (CLPRO409-2) sterile filtered white lyophilized (freeze-dried) powder. The protein was lyophilized from a concentrated solution (1mg/ml) containing no additives.

Solubility: It is recommended to reconstitute the lyophilized S100A6 in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Stability: Lyophilized S100A6 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution S100A6 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).

Please avoid freeze-thaw cycles.

Purity:

Greater than 90.0% as determined by:

Analysis by reducing and non-reducing SDS-PAGE Coomassie.

Applications:

Immunoblots.

Absorption.

Immunohistochemistry

FOR RESEARCH USE ONLY

JV 11/03/06

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com