



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Recombinant Human Bcl-2 (rHuBcl-2)

CLPRO411
CLPRO411-2
Lot:

Description: Recombinant Human Bcl-2 produced in E.Coli is a single, non-glycosylated polypeptide chain containing amino acids 1-218. Recombinant Bcl-2 is expressed as His-Tag fusion protein and purified by proprietary chromatographic techniques.

Source: *E. Coli*

Presentation: 1µg (CLPRO411) or 5µg (CLPRO411-2), sterile filtered colorless solution. The protein contains 10mM Tris-HCL pH 8, 1mM EDTA and 250mM NaCl.

Stability: Recombinant Bcl-2 although stable at 12°C for 1 week, should be stored desiccated below -18°C. Please avoid freeze-thaw cycles.

Purity:

Greater than 95.0% as determined by:

- (a) Analysis by RP-HPLC.
- (b) Analysis by reducing and non-reducing SDS-PAGE Silver Stained gel.

Endotoxin: Less than 0.1 ng/µg (IEU/µg) of Bcl-2.

Applications: Input marker or positive control (Western Blotting).
Function study (binding assay).

Laboratory Reagent For Research Use Only

JV 05/30/07

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com