



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Recombinant Human Recoverin

CLPRO441
CLPRO441-2
CLPRO441-3

Introduction: Recoverin is a heterogeneously acylated calcium-binding and intracellular signal transduction protein in the photoreceptor cells of retina. Recoverin contains four EF-hands, of which two bind Ca. Ca-induced extrusion of the acyl group from a hydrophobic cleft in the protein drives the translocation of recoverin from solution to the disc membrane. Recently, recoverin is a detectable serologic protein that is expressed in patients with cancer-associated retinopathy, a paraneoplastic syndrome.

Description: Recombinant human recoverin produced in *E.Coli* is a single, non-glycosylated polypeptide chain containing 200 amino acids & having a molecular mass of 23 kDa.

Synonyms: RCV1, Cancer-associated retinopathy protein, Protein CAR, RCVRN

Source: *Escherichia Coli*.

Presentation: 10 µg (**CLPRO441**), 50 µg (**CLPRO441-2**), or 1 mg (**CLPRO441-3**), sterile filtered colorless solution. The protein (1mg/ml) contains 20 mM Tris, pH 8.0, 1 mM EDTA, 2 mM MgCl₂ and 10% glycerol.

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles. For long term storage it is recommended to add a carrier protein (0.1% HAS or BSA).

Purity:

Greater than 95.0% as determined by:

- (a) Analysis by RP-HPLC.
- (b) Analysis by SDS-PAGE.

Sequence: MGNSKSGALS KEILEELQLN TKFSEELCS WYQSFLKDCP TGRITQQQFQ SIYAKFFPDT
DPKAYAQHVF RSFDSNLDGT LDFKEYVIAL HMTTAGKTNQ KLEWAFSLYD
VDGNGTISKNEVLEIVMAIF KMITPEDVKL LPDDENTPEK RAEKIWKYFG KNDDDKLTEK
EFIEGTLANK EILRLIQFEP QKVKEKMKNA

Laboratory Reagent For Research Use Only

MW 03/04/15

Visit our website for your local distributor.

CEDARLANE® 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com