



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Alpha Actinin (α -Actinin)

CLPRO518
CLPRO518-2
CLPRO518-3
Lot:

Description: Ultra pure Alpha Actinin having a Molecular mass of 95,000 dalton.

Source: Chicken Gizzard

Presentation: 10 μ g (CLPRO518), 50 μ g (CLPRO518-2), or 1 mg (CLPRO518-3), sterile filtered white lyophilized (freeze-dried) powder. The protein was lyophilized from a 1mg/ml solution containing 10 mM Tris acetate buffer pH 7.6, 0.1 mM EDTA, 2 mM DTT, and 20 mM NaCl.

Solubility: It is recommended to reconstitute the lyophilized α -Actinin in sterile 18M Ω -cm H₂O not less than 100 μ g/ml, which can then be further diluted to other aqueous solutions.

Stability: Lyophilized α -Actinin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution α -Actinin should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). **Please prevent freeze-thaw cycles.**

Applications:

Protein standard in 1D and 2D SDS gel electrophoresis
Immunoassays
Immunization

Purity:

Greater than 95.0% as determined by:
(a) Analysis by RP-HPLC.
(b) Analysis by SDS-PAGE.

Laboratory Reagent For Research Use Only.

JV 11/29/07

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com