



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Human Plasminogen Activator Inhibitor-1 (PAI-1)

CLPRO524

CLPRO524-2

Lot: 606HPAI101

Description: Constitutively active human plasminogen activator inhibitor 1, stable mutant 14-1B1 having a molecular mass of 43 kDa. This human form of PAI-1 contains the following four mutations: K154T, Q139L, M354I and N150H. These mutations combine to confer great stability to the otherwise labile molecule essentially locking it into the active conformation.

Presentation: 2 µg (CLPRO524) or 10 µg (CLPRO524-2), sterile clear solution. The protein solution (2.5 mg/ml) contains 0.05 M sodium phosphate 0.1 M NaCl and 1 mM EDTA (pH 6.6).

Stability: Human Plasminogen Activator Inhibitor-1 should be stored below -20°C. Please avoid freeze-thaw cycles.

Applications:

Protein standard in 1D and 2D SDS gel electrophoresis

Immunoassays

Immunization

ELISA

Binding Studies

Purity:

Greater than 95.0% as determined by:

(a) Analysis by RP-HPLC.

(b) Analysis by SDS-PAGE.

For Research Use Only

JV 11/01/06

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com