



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Technically  
Speaking

**CEDARLANE®**  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Delivering You Today's Innovations for  
the Science of Tomorrow™

## Human High Density Lipoprotein (HDL)

CLPRO559  
CLPRO559-2  
CLPRO559-3  
Lot:

**Introduction:** HDL is a complex of lipids and proteins in approximately equal amounts that functions as a transporter of cholesterol in the blood. HDL is the smallest of the lipoprotein particles; it's the densest because it contains the highest proportion of protein. The liver synthesizes these lipoproteins as complexes of apolipoproteins and phospholipid, which bear a resemblance to cholesterol-free flattened spherical lipoprotein particles. They pick up cholesterol; carry it internally, from cells by interaction with the ATP Binding Cassette Transporter A1. Lecithin-cholesterol acyltransferase (plasma enzyme) converts the free cholesterol into cholesteryl ester (a more hydrophobic form of cholesterol) which is then sequestered into the core of the lipoprotein particle ultimately making the newly synthesized HDL spherical. They increase in size as they pass through the bloodstream and integrate more cholesterol and phospholipid molecules from cells and other lipoproteins, for instance by the interaction with the ABCG1 transporter and the phospholipid transport protein. HDL transports their cholesterol generally to the liver or steroidogenic organs such as adrenals, ovary and testes by direct and indirect pathways. The release of HDL cholesterol to adrenals, ovaries and testes is important for the synthesis of steroid hormones. Since the triglycerides are not stable in HDL, they're degraded by hepatic lipase in order that finally small HDL particles are left which resume the uptake of cholesterol from cells. The cholesterol supplied to the liver is excreted into the bile and consequently intestine either directly or indirectly after conversion into bile acids. High levels of HDL are associated with a decreased risk of atherosclerosis and coronary heart disease.

**Description:** Human High Density Lipoprotein (HDL) produced in Human plasma.

**Synonyms:** High Density Lipoprotein, HDL.

**Source:** Human plasma.

**Presentation:** 10 mg (CLPRO559), 100 mg (CLPRO559-2), or 1 g (CLPRO559-3), sterile filtered clear solution.

**Stability:** Human HDL although stable at 4°C for 1 week, should be stored desiccated below -15°C. (short term i.e. < 3 months) and below -70°C for long term.

**Purity:** Cholesterol level >200 mg/l. <1% LDL.

**Laboratory Reagent For Research Use Only**

JV 04/18/11

Visit our website for your local distributor.

**CEDARLANE®**



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)