



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Recombinant Human Charcot-Leyden Crystal Protein

CLPRO744
CLPRO744-2
CLPRO744-3

Introduction: Eosinophil lysophospholipase (CLC) acts on biological membranes to regulate the multifunctional lysophospholipids. CLC is a lysophospholipase expressed in eosinophils and basophils. CLC hydrolyzes lysophosphatidylcholine to glycerophosphocholine and a free fatty acid. The CLC protein may possess carbohydrate or IgE-binding activities. CLC is both structurally and functionally related to the galectin family of beta-galactoside binding proteins. CLC may be linked with inflammation and some myeloid leukemias.

Description: LGALS10 produced in *E.Coli* is a single, non-glycosylated polypeptide chain containing 162 amino acids (1-142 a.a.) and having a molecular mass of 18.6kDa.
LGALS10 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Synonyms: Eosinophil lysophospholipase, Charcot-Leyden crystal protein, CLC, Galectin-10, Gal-10, Lysolecithin acylhydrolase, GAL10, LGALS10, LGALS10A.

Source: *Escherichia Coli*.

Presentation: 5 µg (CLPRO744), 20 µg (CLPRO744-2), or 1 mg (CLPRO744-3), sterile filtered clear solution. Galectin-10 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 10% glycerol and 0.1M NaCl.

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Purity: Greater than 90.0% as determined by SDS-PAGE.

AA Sequence: MGSSHHHHHH SSGLVPRGSH MSLLPVPYTE AASLSTGSTV TIKGRPLACF LNEPYLQVDF
HTEMKEESDI VFHFQVCFGR RVVMNSREYG AWKQQVESKN MPFQDGQEFE LSISVLPDKY
QVMVNGQSSY TFDHRIKPEA VKMVQVWRDI SLTKFNVSYL KR.

Laboratory Reagent For Research Use Only

JK 02/08/18

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com