



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Human Aryl Hydrocarbon Receptor Interacting Protein

CLPRO886
CLPRO886-2
CLPRO886-3
Lot:

Introduction: Aryl hydrocarbon receptor (AHR) interacting protein (AIP) may have a positive role in AHR-mediated (aromatic hydrocarbon receptor) signaling, probably by influencing its receptivity for ligand and/or its nuclear targeting. AIP is a Cellular negative regulator of the hepatitis B virus (HBV) X protein. Furthermore, AIP is a ubiquitously expressed protein, which binds to HSP 90 and AHR through a highly conserved carboxy-terminal tetratricopeptide repeat domain.

Description: AIP Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 350 amino acids (1-330 a.a.) and having a molecular mass of 39.8kDa. AIP is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Synonyms: AH receptor-interacting protein, AIP, Aryl-hydrocarbon receptor-interacting protein, HBV X-associated protein 2, XAP-2, Immunophilin homolog ARA9, XAP2, ARA9, FKBP16, FKBP37, SMTPHN.

Source: Escherichia Coli.

Presentation: 5 µg (CLPRO886), 20 µg (CLPRO886-2), or 1 mg (CLPRO886-3), sterile filtered clear solution. AIP protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 1mM DTT and 10% glycerol.

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Purity: Greater than 90.0% as determined by SDS-PAGE.

AA Sequence: MGSSHHHHHSSGLVPRGSHMADIARLREDGIQKRVIQEGRGELPDFQDGTKATFHYRT
LHSDDEGTVLDDSRARGKPMELIIGKKFKLPVWETIVCTMREGIEAQFLCDIKHVVLVPLVAKSLRNIAVG
KDPLEGQRHCCGVAQMREHSSLGHADLDALQQNPQPLIFHMEMLKVESPGTYQQDPWAMTDEEKAKAV
PLIHQEGNRLYREGHVKEAAKYYDAIACLNLMQKEQPGSPEWQLDKQITPLLLNYCQCKLVVEEYEE
VLDHCSSILNKYDDNVKAYFKRGKAHAHVWNAQEAQADFAKVLELDPALAPVVSRELRALEARIRQKDE
EDKARFRGIFSH.

Laboratory Reagent For Research Use Only

JK 06/26/14

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com