



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Recombinant Human Ribosomal Protein S3

CLPRO954
CLPRO954-2
CLPRO954-3

Introduction: 40S ribosomal protein S3 is a member of the S3P family of ribosomal proteins. RPS3 facilitates the formation of the mRNA binding site 3' of the codon in the decoding site. Moreover, RPS3 takes part in DNA damage recognition due to its affinity for abasic sites and 7, 8-dihydro-8-oxoguanine residues and its collaboration with human base excision repair proteins OGG1 and Ref-1. Expression of RPS3 is higher in colon adenocarcinomas and adenomatous polyps compared to adjacent normal colonic mucosa.

Description: RPS3 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 263 amino acids (1-243) and having a molecular mass of 28.8 kDa. RPS3 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Synonyms: S3, Ribosomal protein S3, 40S ribosomal protein S3, IMR-90 ribosomal protein S3, MGC87870, FLJ26283, FLJ27450.

Source: *Escherichia Coli*.

Presentation: 2 µg (CLPRO954), 10 µg (CLPRO954-2), or 1 mg (CLPRO954-3), sterile filtered colorless solution. The RPS3 solution (1mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 0.1M NaCl and 10% glycerol.

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Purity: Greater than 90% as determined by SDS-PAGE.

Amino acid sequence: MGSSHHHHHH SGLVPRGSH MAVQISKRRK FVADGIFKAE LNEFLTRELA EDGYSQVEVR VTPTRTEIII LATRTQNVLG EKGRRIRELT AVVQKRFGFP EGSVELYAEK VATRGLCAIA QAESLRYKLL GGLAVRRACY GVLRFIMESG AKGCEVVVSG KLRGQRAKSM KFVDGLMIHS GDPVNYVVD AVRHVLLRQG VLGKVKIML PWDPTGKIGP KKPLPDHVS IVEPKDEILPT TPSEQKGGK PEPPAMPQPV PTA.

Laboratory Reagent For Research Use Only

JK 11/15/16

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com