



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Non-Immune Rabbit IgG Purified IgG

CLRB00

Lot: IG2989R1

Exp: 11/2019

Presentation: **1.0 ml** of purified immunoglobulin from normal (non-immunized) rabbit serum. The IgG concentration is **10 mg/ml** in 10 mM HEPES, pH 7.4, and 0.15 M NaCl, containing 50% glycerol.

Storage: Store at -20°C. Avoid storage in frost-free freezers. Allow product to warm to room temperature and gently mix immediately before use.

Specificity: No reaction was observed with normal human plasma using immunoelectrophoresis and immunodiffusion techniques.

Application: This product can be used as a negative control in procedures where rabbit antibodies are used, such as ELISA and immunoblotting.

Optimum conditions for use of this immunoglobulin in analytical procedures should be established in each laboratory

Laboratory Reagent For Research Use Only

MW 07/22/15

Visit our website for your local distributor.

CEDARLANE® 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com