



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Recombinant Rubella Virus E1 Mosaic

CLRUB291
CLRUB291-2
CLRUB291-3
Lot:

Introduction: Rubella virus is an enveloped positive-strand RNA virus of the family TOGAVIRIDAE: Virions are composed of three structural proteins: a capsid and two membrane-spanning glycoproteins, E2 and E1. During virus assembly, the capsid interacts with genomic RNA to form nucleocapsids. The rubella virus (RV) structural proteins: capsid, E2, and E1 are synthesized as a polyprotein precursor. The signal peptide that initiates translocation of E2 into the lumen of the endoplasmic reticulum remains attached to the carboxy terminus of the capsid protein after cleavage by signal peptidase.

Description: The *E.Coli* derived recombinant protein contains the Rubella Virus E1 regions, 157-176, 374-390, 213-239 amino acids.

Purification Method: Rubella protein was purified by proprietary chromatographic technique.

Purity: Rubella protein is >95% pure as determined by 10% PAGE (Coomassie staining).

Presentation: 100 µg (CLRUB291), 500 µg (CLRUB291-2), or 1000 µg (CLRUB291-3) in solution containing 20 mM imidazol, 8 M urea and 0.3 M NaCl.

Storage: Rubella Protein is shipped at ambient temperature. Upon arrival, Store at -20°C.

Stability: Five years frozen. One month in solution at room temperature.

Specificity: Immunoreactive with sera of Rubella Virus infected individuals.

Applications: Rubella antigen is suitable for ELISA and Western blots, excellent antigen for detection of Rubella Virus with minimal specificity problems.

Laboratory Reagent For Research Use Only

JV 09/30/09

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com