



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Normal Rabbit Serum, Rat Cell Adsorbed, IgG Fraction

CLSDG403-R

Lot: 6507RAD-GG

Description: Rabbit serum is prepared from normal, non-immunized rabbits. It is processed from whole blood. This normal rabbit serum is first adsorbed with rat red blood cells, thymocytes, and spleen cells. The IgG fraction is then made from this adsorbed normal rabbit serum. These are the same procedures used to make our IgG fractions of adsorbed rabbit antisera to rat macrophage, neutrophil, and thrombocyte cells (Cat# CLADG51240, CLADG51140 and CLADG51440 respectively). The IgG fractions of adsorbed normal rabbit serum is the control serum when using the above three IgG fractions of adsorbed antisera in immunological tests in vivo or in vitro.

Presentation: 20 mg (5.0 mg IgG/ml in sterile PBS).

Storage: Store at -20°C.

Shelf-Life: When stored properly, this product will keep its usefulness for 2 years.

Preservatives: None (0.2 µm filtered)

Sterility test: None

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Laboratory Reagent For Research Use Only

JK 01/23/17

Visit our website for your local distributor.

CEDARLANE® 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com