



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Mouse Anti-Human Granulocyte Associated Antigen, CD15 Monoclonal Antibody

CLSG36567

Lot:

Immunogen:	Purified neutrophils from normal human peripheral blood.
Clone:	C3D-1
Isotype:	IgM, Kappa
Format:	This antibody is supplied as tissue culture supernatant containing sodium azide as preservative.
Titer/Working Dilution:	This antibody may be diluted to a titer of 1:10-1:25 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.
Staining Protocol:	We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. <u>Suitable for formalin/paraffin tissues.</u>
Specificity:	The specificity of this antibody is similar to antibodies Leu-M1, Tü 9 and VIM-D5, recognizing an epitope involving the carbohydrate sequence 3-fucosyl-N-acetylacetosamine. This antibody reacts with an antigen located on mature granulocytes and on Hodgkin and Reed-Sternberg cells.
Positive Control:	Tonsil
Cellular Localization:	Cell membrane
Storage:	Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.
References:	i) Stein et al. J Clin Pathol 86: 292, 1986. ii) Stein et al. Int J Cancer 29: 283, 1982. iii) Stein et al. Int J Cancer 30: 445, 1982. iv) Dorfman et al. Am J Pathol 123: 508, 1986.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

JV 12/17/05

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com