



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Technically
Speaking

CEDARLANE®
www.cedarlanelabs.com



Delivering You Today's Innovations for
the Science of Tomorrow™

Mouse Anti-Human B Cell, CD45RA Monoclonal Antibody

CLSG36597
CLSG36597-05
Lot:

Immunogen: BALB/C mice were immunized with hairy cell leukemia cells.

Clone: 4KB5

Isotype: IgG1, Kappa

Format: 1 ml (CLSG36597), or 0.5 ml (CLSG36597-05). This antibody is supplied as tissue culture supernatant containing sodium azide as a preservative.

Titer/Working Dilution: This antibody may be diluted to a titre of 1:10-1:25 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. Suitable for formalin fixed, paraffin embedded tissue sections.

Specificity: This antibody reacts with most B cells in tissue sections and peripheral blood as well as small population of T cells and onocytes. Clone 4KB5 appears similar to two other anti-CD45R antibodies.

Positive Control: Tonsil

Cellular Localization: Cell membrane

Storage: Store at 2-8°C. Do not use beyond the expiration date stated on the label.

References:
i) Dalchau et al. J Exp Med 153: 753, 1981.
ii) Poppema et al. Am J Pathol 127: 418, 1987.

Laboratory Reagent For Research Use Only

JV 04/12/11

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com