



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Mouse Anti-Human Cytokeratin 10

CLSG36674

Lot:

Description:

This monoclonal antibody recognizes keratin polypeptide of 56.5 kD. It stains suprabasal cells in the epidermis and thymic Hassall's bodies. Simple and glandular epithelia, transitional epithelium and adenocarcinomas are negative with this antibody.

Presentation:

1.0 ml, purified immunoglobulin containing 0.05% sodium azide (NaN_3) as a preservative.

Immunogen:

BALB/C mice were immunized with detergent insoluble extract of human epidermis.

Isotype:

Mouse IgG₁ kappa

Clone:

DE-K10

Application:

This antibody is suitable for use in immunohistochemistry. We recommend a starting dilution of 1:25 – 1:50 with an ABC method, and an incubation period of 30 minutes at room temperature. Enzymatic predigestion of paraffin embedded tissues is **not** recommended prior to the immunostaining, as it may reduce the staining. This antibody is also suitable in Western Blotting (1:50 to 1:100).

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

phone: (905) 878-8891 • fax: (905) 878-7800

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

For optimal results in various applications, it is recommended that each investigator determine dilutions and incubation period appropriate for individual use.

Positive Control: Skin

Cellular Localization: Cytoplasmic

Storage: 4°C. Check label for expiry date.

References:

1. Ivanyi *et al.* J. Pathol. 159: 7, 1989.
2. Ivanyi *et al.* Differentiation 42: 124, 1989.
3. Ivanyi *et al.* Cancer Res. 50: 5143, 1990.

FOR RESEARCH USE ONLY

TN 05/13/05