



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Mouse Anti-Human Cytokeratin 16

CLSG36680

Lot:

Description:

This monoclonal antibody is specific to a 48 kD protein, which is identified as cytokeratin 16. Cytokeratin 16 is expressed in keratinocytes, which are undergoing rapid turnover in the suprabasal region. High concordance was found between the carcinomas immunostaining with the basal cell and the hyper proliferation related keratins.

Presentation:

1.0 ml, purified immunoglobulin fraction containing 0.05% sodium azide (NaN_3) as a preservative.

Immunogen:

BALB/C mice were injected with a synthetic peptide from the C-terminus of human cytokeratin 16.

Isotype:

Mouse IgG₁

Clone:

LL025

Application:

This antibody is suitable for use in immunohistochemistry. We recommend a starting dilution of 1:25 – 1:75 with an ABC method, and an incubation period of 30 minutes at room temperature. Formalin-fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10mM citrate buffer, pH 6.0 prior to immunostaining.

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

For optimal results in various applications, it is recommended that each investigator determine dilutions and incubation period appropriate for individual use.

Positive Control: Skin

Cellular Localization: Cytoplasmic

Storage: 4°C. Check label for expiry date.

References:

1. Markey *et al.* Br. J. Dermatol. 126: 154, 1992
2. Wetzels *et al.* Am. J. Pathol. 138:751, 1991

FOR RESEARCH USE ONLY

TN 05/12/05