



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Mouse Anti-Skeletal Myosin (Fast) Monoclonal Antibody

CLSG36794
CLSG36794-05
Lot:

Immunogen: Rabbit muscle myosin.

Clone: MY32

Isotype: IgG1

Format: 1 ml (**CLSG36794**) or 0.5 ml (**CLSG36794-05**). This antibody is supplied as diluted ascites containing sodium azide as a preservative.

Titer/Working Dilution: This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed.

Specificity: This antibody recognizes the heavy chain of human skeletal myosin. It does not stain cardiac or smooth muscle myosin. It stains the fast twitch (type II) isomyosin molecule and cross reacts with myosin from rabbit, rat, mouse, bovin, chicken, and guinea pig.

Positive Control: Skeletal muscle

Cellular Localization: Cytoplasmic

Storage: Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

Laboratory Reagent For Research Use Only

BA 03/27/12

Visit our website for your local distributor.

CEDARLANE® 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com