



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Place your order with CEDARLANE® or your local distributor.

*Please contact CEDARLANE® for lot specific information.*

### **Mouse Anti-Human p300/CBP**

**CLSG36815**

**Lot:**

**Description:** This monoclonal antibody is directed against human p300/CBP, and detects a 300 kD protein. P300 and CBP (CREB-binding protein) are homologous to each other. P300 is a transcriptional adapter protein, which is characterized by three cysteine and histidine-rich regions. Protein kinase A mediated CREB phosphorylation results in the binding of CREB to a 265 kD nuclear protein designated CBP. This antibody cross reacts with monkey, mink, mouse and rat.

**Presentation:** 1 ml, purified ascites containing sodium azide ( $\text{NaN}_3$ ) as a preservative.

**Immunogen:** BALB/C mice were immunized with E1A affinity purified p300 protein.

**Isotype:** Mouse IgG<sub>2b</sub> kappa

**Clone:** NM11

**Applications:** This antibody is suitable for use in immunohistochemistry. We recommend a starting dilution of 1:25 – 1:50 with an ABC method, and an incubation period of 30 minutes at room temperature. Formalin-fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10mM citrate buffer, pH 6.0 prior to immunostaining.

*Continued Overleaf...*

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
**in North America**

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)

For optimal results in various applications, it is recommended that each investigator determine dilutions and incubation period appropriate for individual use.

**Positive Control:** Tonsil

**Cellular Localization:** Nuclear

**Storage:** 4°C. Check label for expiry date.

**References:**

1. Dallas *et al.* J. Virol. 71: 1726, 1997.
2. Dallas *et al.* Hybridoma 16:273, 1997.

**FOR RESEARCH USE ONLY**

TN 03/24/05