



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.
Please contact CEDARLANE® for lot specific information.

Purified Rabbit Anti-Calretinin Polyclonal Antibody

CLSG36897

Lot:

Immunogen: Recombinant full length mouse Calretinin protein.

Host: Rabbit

Format: 1 ml of affinity purified rabbit antiserum against Calretinin containing sodium azide as preservative.

Titer/Working Dilution: This antibody may be diluted to a titer of 1:50-1:100 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature or overnight at 4 degrees Celsius. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining.

Specificity: The antibody is specific against a calcium binding protein called calretinin. This antibody reacts specifically with calretinin in tissue originating from human, rat and mouse. Calretinin is an intracellular calcium-binding protein belonging to the troponin C superfamily. Calretinin has been detected in mesotheliomas and some pulmonary adenocarcinomas.

Positive Control: Mesothelioma

Cellular Localization: Cytoplasmic

Storage: Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.

Continued Overleaf...

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

References:

Dogliani et al. Am J Surg Pathol 20: 1037, 1996.
Rogers J. J Cell Biol 105: 1343, 1987.
Schwaller et al. Cell Calcium 14: 639, 1993.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

JV 05/03/06