



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically
Speaking*

CEDARLANE® 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Purified Rabbit anti-Human IgA Polyclonal Antibody

CLSG36948
CLSG36948-05
Lot: I507

Immunogen:	IgA isolated from a pool of normal human plasma.
Host:	Rabbit
Format:	1 ml (CLSG36948) or 0.5 ml (CLSG36948-05) affinity purified immunoglobulin fraction of rabbit antiserum against human IgA containing sodium azide as preservative.
Titer/Working Dilution:	This antibody may be diluted to a titer of 1:25-1:100 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.
Staining Protocol:	This antibody is to be used with FFPE tissue sections pre-treated in citrate buffer, pH 6.0. We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user.
Specificity:	This antibody reacts with the alpha chains of human IgA. This antibody has been adsorbed to a very high level of specificity.
Positive Control:	Tonsil
Cellular Localization:	Cytoplasmic
Storage:	Store the antibody at 2-8°C.

For Research Use Only. Not for Diagnostic or Therapeutic Work

MW 09/11/18

Visit our website for your local distributor.

CEDARLANE® 
www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com