



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Purified Rabbit anti-Human Ki-67 Antigen Polyclonal Antibody

CLSG36954

Lot:

Immunogen: Synthetic peptide from 62 base pair region of the human Ki-67 antigen.

Host: Rabbit

Format: 1 ml of purified immunoglobulin fraction of rabbit antiserum against human Ki-67 antigen containing sodium azide as preservative.

Titer/Working Dilution: This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining.

Specificity: This antibody reacts with a nuclear antigen present in proliferating human cells. This antibody cross-reacts with a range of mammalian species. It can be used to evaluate the Ki-67 labeling index in various tumors.

Positive Control: Tonsil

Cellular Localization: Nuclear

Storage: Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.

References: i) Gerdes et al. Am J Pathol 138: 867, 1991.

ii) Wintzer et al. Cancer 67: 421, 1991.

iii) Brown et al. Histopathology 17: 489, 1990

iv) Gerdes J. Seminars in Cancer Biol 1: 199, 1990.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

JV 07/19/06

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com