



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.
Please contact CEDARLANE® for lot specific information.

Rabbit Anti-Human Kappa Light Chain FITC-Conjugated

CLSG80087
Lot: G858

Immunogen:	Pool of human kappa light chain isolated from the urine of patients with Bence Jones Proteinuria.
Format:	1 ml purified immunoglobulin fraction of rabbit antiserum conjugated with fluorescein isothiocyanate isomer 1, containing sodium azide as preservative.
Titer/Working Dilution:	This antibody may be diluted to a titer of 1:20-1:40 for immunofluorescence technique. However, depending on the staining conditions employed, we suggest that the final dilution should be determined by the user.
Staining Protocol:	We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions optimal incubation conditions should be determined by the end user.
Specificity:	Traces of contaminating antibodies have been removed by solid phase absorption. Only the kappa precipitation line appeared in cross immunoelectrophoresis. This antibody is well suited for tissue section staining but may be used for other applications also.
Storage:	Store the antibody at 2-8°C in dark. Do not use it beyond the expiration date stated on the label.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

JV 02/14/08

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com