



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

Technically  
Speaking

**CEDARLANE**<sup>®</sup>  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Delivering You Today's Innovations for  
the Science of Tomorrow™

## FITC Rabbit Anti-Human C3c Complement

**CLSG80089**

**Lot:**

**Immunogen:**

C3c Complement isolated from complement-activated human serum.

**Format:**

1ml purified immunoglobulin fraction of rabbit antiserum conjugated with fluorescein isothiocyanate isomer 1, containing sodium azide as a preservative.

**Titer/Working Dilution:**

This antibody may be diluted to a titer of 1:75-1:100 for immunofluorescence technique. The final dilution should be determined by the user based upon the staining conditions employed.

**Staining Protocol:**

We suggest an incubation period of 30-60 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed.

**Specificity:**

This antibody reacts with the human C3c complement and the C3c part of C3 and C3b. It cross reacts with the C3c complement from cat, cow, dog, goat, guinea pig, kangaroo, mink, mouse, rat, sheep, and swine. Traces of contaminating antibodies have been removed by solid phase absorption with human plasma proteins. This antibody is well suited for tissue section staining but may be also used for other applications.

**Storage:**

Antibody is light sensitive. Store at 2-8°C in the dark. Do not use beyond the expiration date stated on the label.

**Laboratory Reagent For Research Use Only**

JV 11/25/10

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup>



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)