



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.  
Please contact CEDARLANE® for lot specific information.

## Rabbit Anti-Human C4c Complement FITC-Conjugated

**CLSG80091**

**Lot:** H370

**Immunogen:**

C4c Complement isolated from complement-activated human serum.

**Format:**

1ml of purified immunoglobulin fraction of rabbit antiserum conjugated with fluorescein isothiocyanate isomer 1, containing sodium azide (NaN<sub>3</sub>) as preservative.

**Titer/Working Dilution:**

This antibody may be diluted to a titer of 1:50-1:100 for immunofluorescence technique. However, depending on the staining conditions employed, we suggest that the final dilution should be determined by the user.

**Staining Protocol:**

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions optimal incubation conditions should be determined by the end user.

**Specificity:**

This antibody reacts with human C4, C4b and C4c but does not react with the C4d epitope. Traces of contaminating antibodies have been removed by solid phase absorption with human plasma proteins. This antibody is well suited for tissue section staining but may be used for other applications also.

**Storage:**

Store the antibody at 2-8°C in dark. Do not use it beyond the expiration date stated on the label.

**For Research Use Only. Not for Diagnostic or Therapeutic Work.**

JV 01/21/07

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
**in North America**

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)