



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Delivering You Today's Innovations for  
the Science of Tomorrow<sup>™</sup>

## Rabbit Anti-Human Fibrinogen FITC-Conjugated

**CLSG80092**

**Lot:** L306

**Immunogen:**

Fibrinogen isolated from human plasma.

**Format:**

**1ml** of purified immunoglobulin fraction of rabbit antiserum conjugated with fluorescein isothiocyanate isomer 1, containing sodium azide (NaN<sub>3</sub>) as preservative.

**Titer/Working Dilution:**

This antibody may be diluted to a titer of 1:20-1:40 for immunofluorescence technique. However, depending on the staining conditions employed, we suggest that the final dilution should be determined by the user.

**Staining Protocol:**

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions optimal incubation conditions should be determined by the end user.

**Specificity:**

This antibody reacts with native fibrinogen as well as with the fibrinogen fragments D, E, X, and Y. Traces of contaminating antibodies have been removed by solid phase absorption with human plasma proteins. This antibody is well suited for tissue section staining but may be used for other applications also.

**Storage:**

Store the antibody at 4-8°C in dark. Do not use it beyond the expiration date stated on the label.

**Laboratory Reagent For Research Use Only**

JV 05/19/09

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup>



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)