



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE®
www.cedarlanelabs.com



Delivering You Today's Innovations for
the Science of Tomorrow™

Mouse Anti-Pan Muscle Actin Monoclonal Antibody

CLSG80217
CLSG80217-05
Lot:

Immunogen: BALB/C mice were injected with N-terminal decapeptide of alpha-smooth muscle actin for 1A4 and with purified rabbit striated muscle actin for 5C5.

Clone: 1A4 & 5C5

Isotype: IgG2a, kappa & IgM

Format: 1 ml (**CLSG80217**) or 0.5 ml (**CLSG80217-05**). This antibody is supplied as diluted ascites containing sodium azide as a preservative.

Titer/Working Dilution: This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. Suitable for formalin fixed, paraffin embedded tissue sections.

Specificity: This pan actin is useful in identifying tumors arising from muscle cells, i.e. leiomyosarcoma and rhabdomyosarcoma. It cross-reacts with human, sheep, bovine, rabbit, guinea pig, rat, frog, snake, and carp.

Positive Control: Skeletal muscle

Cellular Localization: Cytoplasmic

Storage: Store at 2-8°C. Do not use beyond the expiration date stated on the label.

References:
i) Babai et al. Virchow Arch 55: 263, 1988.
ii) Schurch et al. Amer J Pathol 128:91, 1987.
iii) Skalli et al. Am Amer J Pathol 130: 515, 1988.

Laboratory Reagent For Research Use Only

JV 05/19/11

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com