



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Mouse Anti-Human Alpha-Smooth Muscle Actin Monoclonal Antibody

CLSG80258

Lot#:

Immunogen:

BALB/C mice were injected with N-terminal decapeptide of alpha-smooth muscle actin.

Size:

6 ml

Clone:

1A4

Isotype:

IgG2a, Kappa

Format:

This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.

Staining Protocol:

We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Suitable for formalin/paraffin tissues.

Specificity:

This antibody is specific to alpha-smooth isoform of actin. It reacts with smooth muscle cells of vessels and different parenchymes. This antibody does not cross-react with beta and gamma-cytoplasmic, alpha-sarcomeric and alpha-myocardial actin isoforms.

Positive Control:

Leiomyoma

Cellular Localization:

Cytoplasmic

Storage:

Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.

References:

- i) Skalli et al. Am J Cell Biol 103: 2787, 1986.
- ii) Skalli et al. Lab Invest 60: 275, 1989.
- iii) Skalli et al. J Histo Chem 37: 315, 1989.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

TN 4/19/2006

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com