



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE®
www.cedarlanelabs.com



Delivering You Today's Innovations for
the Science of Tomorrow™

Mouse anti-Human Progesterone Receptor Pre-Diluted Monoclonal Antibody

CLSG80278
CLSG80278-1
Lot:

- Immunogen:** Recombinant human progesterone receptor protein.
- Clone:** PGR-1A6
- Isotype:** IgG1, Kappa
- Format:** 5 ml (CLSG80278), or 6 ml (CLSG80278-1). This antibody is supplied in a ready-to-use format.
- Specificity:** This antibody is specific to the human progesterone receptor protein. Progesterone receptor is an estrogen regulated protein and its expression indicates a responsive estrogen receptor pathway.
- Positive Control:** Breast carcinoma
- Cellular Localization:** Nuclear
- Storage:** Store at 2-8°C. Do not use beyond the expiration date stated on the label.
- References:**
- i) Mason et al. Cancer Res 43: 2985, 1983.
 - ii) Clark et al. The New England J Med 309: 1343, 1983.
 - iii) Deamant et al. Applied Immunohistochemistry 1: 188, 1993.

Laboratory Reagent For Research Use Only

JV 09/07/10

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com