



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE®
www.cedarlanelabs.com



Delivering You Today's Innovations for
the Science of Tomorrow™

Mouse Anti-Human Cytokeratin (H.M.Weight) Prediluted Monoclonal Antibody

CLSG80321
Lot

Immunogen: BALB/c mice were immunized with solubilized keratin extract from human stratum corneum.

Clone: 348E12

Isotype: IgG1, kappa

Format: **6 ml.** This antibody has been pretitered and quality controlled to work on formalin-fixed, paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Proteolytic treatment is required prior to immunostaining of formalin fixed paraffin embedded tissue sections.

Specificity: This antibody recognizes keratin polypeptide of 68, 58, 56.5 and 50 kDa in extract of stratum corneum. The antibody reacts with squamous, ductal and other complex epithelia. It stains with adenocarcinomas, breast, pancreas, bile duct, salivary gland and transitional cell carcinomas.

Positive Control: Skin

Cellular Localization: Cytoplasmic

Storage: Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

References:
i) Gown et al. J Cell Biol 95: 414, 1982.
ii) Gown et al. Am J Clin Pathol 84: 413, 1985.
iii) O'Malley et al. Virch Arch A Pathol Anat 417: 191, 1990.

Laboratory Reagent For Research Use Only

JV 06/26/09

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com