



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

## Mouse Anti-Human Glial Fibrillary Acidic Protein (GFAP) Monoclonal Antibody

CLSG80325

Lot:

**Immunogen:** BALB/c mice were immunized with glial fibrillary acidic protein from human brain.

**Clone:** 6F2

**Isotype:** IgG1, kappa

**Format:** 6 ml. This antibody has been pre-titered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.

**Staining Protocol:** We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Suitable for formalin/ paraffin tissues.

**Specificity:** This antibody is specific to a 52 kDa intermediate filament protein in brain and spinal cord. It stains astrocytes and CNS ependymal cells but does not stain oligodendrocytes or neurons. This antibody does not cross react with any other intermediate filament proteins. However, it cross reacts with GFAP from mouse and rabbit.

**Positive Control:** Brain

**Cellular Localization:** Cytoplasmic

**Storage:** Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

**References:**

- i) Debus et al. Differentiation 25: 193, 1983.
- ii) Miettinen et al. Ultrastruct Pathol 7: 83, 1984.
- iii) Royds et al. Acta Neuropathol 70: 320, 1986.

**Laboratory Reagent For Research Use Only**

JV 05/30/08

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
**in North America**

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)