



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Mouse Anti-Human CD10 antigen Pre-diluted Monoclonal Antibody

CLSG80351

Lot:

Immunogen: Prokaryotic recombinant protein corresponding to the external domain of the CD10 glycoprotein.

Clone: 56C6

Isotype: IgG1

Format: 6 ml. This antibody has been pre-titered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.

Staining Protocol: We suggest an incubation period of 30 to 60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM Citrate buffer, pH 6.0 prior to immunostaining.

Specificity: This antibody is specific to human CD10 antigen of 100 kDa, also known as common acute lymphocytic leukemia antigen (CALLA). CD10 antigen has been identified on the surface of early lymphoid progenitor cells, immature B cells with in adult bone marrow and germinal center B cells with in lymphoid tissue.

Positive Control: Tonsil

Cellular Localization: Cell membrane

Storage: Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

References:

- i) Carrel et al. Melanoma Res 3: 319, 1993.
- ii) Haralambidou et al. J Clin Pathol 40: 490, 1987.
- iii) Kiyokawa et al. Clin Exp Immunol 79: 322, 1990.
- iv) Scheuermann et al. Leukemia and Lymphoma 18: 385, 1995.

Laboratory Reagent For Research Use Only

JV 11/13/07

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com