



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Rabbit anti-Human Alpha-1-Antichymotrypsin Prediluted Polyclonal Antibody

CLSG80395

Lot:

Immunogen:

Alpha-1-antichymotrypsin isolated from human serum.

Host:

Rabbit

Format:

1 ml. This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.

Staining Protocol:

We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user.
Suitable for formalin/paraffin tissues.

Specificity:

The antibody is specific against human alpha-1-antichymotrypsin. This antibody stains monocytes and macrophages in human tissues and does not stain any other cell type. This antibody shows cross reactivity against wide range of mammalian species.

Positive Control:

Tonsil

Cellular Localization:

Cytoplasmic

Storage:

Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.

References:

1. Meis et al. Am J Clin Pathol 86: 591, 1986.
2. Nathrath et al. Acta Histochem 25: 69, 1982.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

JV 06/19/06

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com